

CMC Considerations for Oncolytic Viral Products

Bo Liang, Ph.D.

CMC Reviewer, Gene Therapy Branch
Division of Cellular and Gene Therapies
Office of Tissues and Advanced Therapies
Center for Biologics Evaluation and Research
Food and Drug Administration

SBIA-Oncology Therapy Development Workshop
March 30-31, 2021


Outline

- Background
- Product manufacturing and testing
- Monitoring shedding of oncolytic virus during clinical trials
- Summary

Oncolytic Virus (OV)

- Infects and lyses cancer cells but not normal cells
- MOAs: Direct tumor lysis and stimulation of immune responses
- Types
 - Naturally occurring
 - Genetically engineered
 - Expression of transgene to enhance the therapeutic effect
 - Modifications that reduce virulence and promote tumor-specific targeting

Examples of OVs



Clinical Landscape of OV_s

- Tumor types: Solid tumors and hematological tumors
- ROAs: Intratumoral, intravenous, intraperitoneal, intravesical, etc.
- Combination therapies: Chemotherapy, radiation, pro-drugs, immune checkpoint inhibitors, etc.





Figure 1 Pie chart showing the distribution of oncolytic viruses by clinical stage. The majority of studies were phase I (n=49; 51%) or not specified (n=29; 30%). There were 6 (6%) phase I/II trials, 11 (11%) phase 2 and only two phase 3 clinical trials.


Product Characterization Increases with Phase of Clinical Investigation




Ensuring Product Safety and Quality



- Qualification of components
- A well-defined process and process controls
- Appropriate product testing



Components




Master and Working Cell Banks

- Origin, derivation history, culture conditions
- Qualification Testing
 - Identity *
 - Sterility (bacteriostasis/fungistasis) *
 - Mycoplasma *
 - Adventitious viral agent testing (AVA): *in vitro* * and *in vivo*
 - Detection of viral particles by Transmission Electron Microscopy (TEM) analysis
 - Retroviruses (reverse transcriptase activity)
 - Bovine and porcine, and other species-specific viruses
 - Viability *
 - Stability
 - Tumorigenicity (required prior to Phase 3 clinical trials)

Master and Working Virus Banks

- Origin, derivation history, culture conditions
- Testing and characterization
 - Sterility (bacteriostasis/fungistasis)
 - Mycoplasma
 - Adventitious viral agent testing: in vivo and in vitro
 - Species-specific pathogens
 - Identity (e.g. by sequencing)
 - Activity
 - Titer/infectivity
 - Differential infection/cytotoxicity
 - Stability
 - Additional testing may be required on a case-by-case basis, e.g., drug sensitivity, transgene expression, etc.




Reagents


- Media
- Serum
- Media supplements (e.g., transferrin, insulin, trypsin, etc.)
- Ancillary materials (e.g. buffer, chemicals, WFI, etc.)
- Reagents used in downstream purification (e.g., nucleases, column materials)

Selection of the highest quality reagents available, establishing vendor/reagent qualification, documentation of the source and derivation is important

Process



Testing



In-Process Testing

- Examples:
 - Cell density
 - Cell viability
 - Bioburden
 - Process/product related impurities
 - Viral titers
- Mycoplasma and *in vitro* AVA testing on Unclarified Bulk Harvest

In vitro AVA Testing

- Challenging due to amplification of the viral product causing cytopathic effects (CPE) in testing cells
- Alternative testing approaches:
 - Use antibodies to neutralize the virus prior to testing (dilution as little as possible)
 - Test on uninfected production control cells (parallel culture)
 - Use advanced and sensitive method, e.g., NGS, as a supplementary approach to the standard AVA testing

Product Testing



- Sterility (CFR 610.12)
- Identity (CFR 610.14)
 - Nucleic acid: e.g., PCR, Sequencing
 - Protein: e.g., Western Blot, ELISA
- Purity (CFR 610.13)
 - Endotoxin
 - Residual contaminants (DNA, protein, culture reagents, column materials, etc.)
 - Particle to Infectivity (PI) ratio
 - Sequence variants
- Potency (CFR 610.10)
 - Direct or indirect measure of biological activity
 - Titer (TCID50, pfu, or viral particle)
 - Transgene expression/function
 - Tumor-specific cell killing
 - Quantitative
- Appearance
- Reference standards

Product Stability



- Long-term storage
 - Appearance
 - Potency (titer)
 - Sterility
 - Quality (degradation products)
- Accelerated & Stress conditions
- Short-term storage: product formulation at clinic and administration (device compatibility)
 - Potency (titer)
- Shipping
 - Potency (titer)
 - Sterility
 - Quality (degradation products)

Shedding

- “Shedding” – release of oncolytic products from the subject through excreta (feces), secretions (urine, saliva, nasopharyngeal fluids, etc.), or the skin (pustules, sores, wounds)
- Shedding is distinct from biodistribution.
- Assess the potential for environmental release and transmission

Shedding

- Product-specific factors should influence the study design
 - Characteristics of the virus (replication competence, persistence/latency, immunogenicity, tropism, etc.)
 - Route of administration, dose, regimen
- Preclinical shedding data may inform the design of human shedding studies
- Collect shedding data in Phase 1 and continue through Phase 2 and Phase 3

Considerations for Shedding Study Design

- Types of samples collected
- Frequency and duration of sample collection
 - Frequent sampling immediately following administration
 - Continue sampling until 3 consecutive data points \leq LOD
 - Sampling to capture 2nd peak of shedding may be necessary
- Methods of detection
 - Detection of product-specific nucleic acids
 - Analysis of positive samples for infectivity/growth

Summary

- Product safety is a key consideration in early phase studies due to the replication competence of the OV.
- Because OVs are complex biologics manufactured using materials of biological origin and terminal sterilization is not possible, product safety should be established through aseptic processing, material qualification, well-defined manufacturing process, and adequate testing.
- The collection of shedding data in preclinical and clinical studies with well-designed study plan for sample collection and the use of qualified methods is an important aspect of clinical development of OV products.

Guidance



- Guidance for Industry: Chemistry, Manufacturing, and Control (CMC) Information for Human Gene Therapy Investigational New Drug Applications (INDs) (January 2020)
- Guidance for Industry: Characterization and Qualification of Cell Substrates and Other Biological Materials Used in the Production of Viral Vaccines for Infectious Disease Indications (February 2010)
- Guidance for Industry: Potency Tests for Cellular and Gene Therapy Products (January 2011)
- Guidance for Industry: Design and Analysis of Shedding Studies for Virus or Bacteria-Based Gene Therapy and Oncolytic Products (August 2015)
- Guidance for Industry: Determining the Need for and Content of Environmental Assessments for Gene Therapies, Vectored Vaccines, and Related Recombinant Viral or Microbial Products (March 2015)
- ICH Considerations on Oncolytic Viruses (September 2009) and General Principles to Address Virus and Vector Shedding (June 2009)

Contact Information

- **Bo Liang**

Gene Therapy Branch,
Division of Cellular & Gene Therapies
OTAT/CBER/FDA
WO Bldg. 71 Room 6018
10903 New Hampshire Ave.
Silver Spring, MD 20993
bo.liang@fda.hhs.gov



- **Regulatory Questions:**

OTAT Main Line – 240 402 8190

Email: OTATRPMS@fda.hhs.gov and Lori.Tull@fda.hhs.gov

- **OTAT Learn Webinar Series:** <http://www.fda.gov/BiologicsBloodVaccines/NewsEvents/ucm232821.htm>

- **CBER website:** www.fda.gov/BiologicsBloodVaccines/default.htm

Phone: 1-800-835-4709 or 240-402-8010

- **Consumer Affairs Branch:** ocod@fda.hhs.gov

- **Manufacturers Assistance and Technical Training Branch:** industry.biologics@fda.gov

- **Follow us on Twitter:** <https://www.twitter.com/fdacber>

